ELECTRON DIFFRACTION EXPERIMENT

Aim: To measure the de Broglie wavelength of electron beam and to compare it's value with the theoretical value obtained from the accelerating potential.

Apparatus: Electron diffraction setup

Theory:

In 1924 deBroglie theorized that not only light posesses both wave and particle properties, but rather particles with mass - such as electrons - do as well. The wavelength of these 'matter waves' - also known as the de Broglie wavelength - can be calculated from Planks constant h divided by the momentum p of the particle.

In the case of electrons it is $\lambda = h/p = h/mv$. The acceleration of electrons in an electron beam gun with the acceleration voltage V_a results in the corresponding de Broglie wavelength.

We have KE = $mv^2/2$ = eV_a , which results to the velocity v = $(2.e.V_a/m)^{1/2}$. Substituting this we get

$$\begin{split} &\lambda = h/p = h/mv = h/[m.(2.e.V_a/m)^{1/2}] = h/(2meV_a)^{1/2} \\ &\lambda = h/(2meV_a)^{1/2} \ldots \ldots (1) \end{split}$$

Tube geometry

This de Broglie hypothesis can be verified experimentally by electron diffraction experiment. This wavelength can be obtained from the tube geometry as follows; From figure we get,

tan
$$(2\theta) = r/(l1 + l2)$$


where l_1 =L-R and l_2 = $(R^2-r^2)^{1/2}$
So $tan(2\theta) = r/(L-R+(R^2-r^2)^{1/2})$
So, $\theta = (\frac{1}{2}). tan^{(-1)}[r/(L-R+(R^2-r^2)^{1/2})].....(2)$

From diffraction condition,

We have for the first diffraction ring,

$$\lambda = d.\sin(\theta)....(3)$$

where d is lattice spacing of graphite and θ is angle of diffraction which is given by equation(2) above. This angle can be measured experimentally.



Procedure for virtual experiment:

- 1) Go to the website https://virtuelle-experimente.de/en/elektronenbeugung/wellenlaenge/quantitativ.php
- 2) Change the slider of the potential and set the required potential.
- 3) Change the slider of the radius and measure the radius of first order ring.
- 4) For calculations use the standard values of Planck's constant, mass of electron and charge of electron, use the value of lattice constant $d = 2.13 \times 10^{-10} \text{m}$ and the respective tube constants as Radius R = 6.35cm, Length L = 12.7cm.
- 5) Using the radius of the ring r, calculate the experimental wavelength using eqn(2) and eqn(3).
- 6) For the given potential, calculate the theoretical wavelength from eqn(1).

Observation table:

| Sr | Accelerating potential | λ (theoretical) | Radius of the | λ (experimental) |
|----|------------------------|-----------------|------------------|------------------|
| | V_a | using eqn(1) | first order ring | using eqn(2&3) |
| | | | r | |
| 1 | 4.5kV | | | |
| 2 | 6kV | | | |
| 3 | 8kV | | | |
| 4 | 10kV | | | |
| 5 | 12kV | | | |

Conclusion:

The theoretical and experimental values of de Broglie wavelength are in good agreement. With increase in the accelerating potential, the de Broglie wavelength decreases as seen from the observation table.